

ART-sRNA CLONING IN B/c VECTORS

Updated on 12/01/2021 by Alberto Carbonell.

Notes:

-Available B/c vectors are listed in Table I at the end of the section.

-AtMIR390-B/c- and AtTAS1c-B/c-based vectors must be propagated in a *ccdB* resistant *E. coli* strain such as DB3.1.

-Alternatively, *BsaI* digestion of the B/c vector and subsequent ligation of the *amiRNA* oligonucleotide insert can be done in separate reactions

1. OLIGONUCLEOTIDE ANNEALING

-Dilute sense oligonucleotide and antisense oligonucleotide in sterile H₂O to a final concentration of 100 μ M.

-If necessary, prepare 1 ml of Oligo Annealing Buffer by mixing the following

<u>Stock</u>	<u>Volume to add</u>	<u>Final Conc.</u>
H ₂ O	613.4 μ L	-
1 M Tris-HCl (pH 7.5)	60 μ L	60 mM Tris-HCl pH7.5
3 M NaCl	166.6 μ L	500 mM NaCl
1 M MgCl ₂	60 μ L	60 mM MgCl ₂
0.1 M DTT	100 μ L	10 mM DTT

Note: Prepare 1 ml aliquots of Oligo Annealing Buffer and store at -20°C.

-Assemble the annealing reaction in a PCR tube as described below:

Forward oligonucleotide (100 μ M)	2 μ L
Reverse oligonucleotide (100 μ M)	2 μ L
<u>Oligo Annealing Buffer</u>	<u>46 μL</u>
Total volume	50 μ L

The final concentration of each oligonucleotide is 4 μ M.

-Use a thermocycler to heat the annealing reaction 5 min at 94°C and then cool down (0.05°C/sec) to 20°C.

-Dilute the annealed oligonucleotides just prior to assembling the digestion-ligation reaction as described below:

Annealed oligonucleotides	3 μ L
<u>dH₂O</u>	<u>37 μL</u>
Total volume	40 μ L

The final concentration of each oligonucleotide is 0.15 μ M.

Note: Do not store the diluted oligonucleotides.

2. DIGESTION-LIGATION REACTION

-Assemble the digestion-ligation reaction as described below:

B/c vector (x ug/uL)	Y μ L (50 ng)
Diluted annealed oligonucleotides	1 μ L
10x T4 DNA ligase buffer	1 μ L
T4 DNA ligase (400 U/ μ L)	1 μ L
<i>Bsa</i> I (10U/ μ L, NEB)	1 μ L
<u>dH₂O</u>	<u>to 10 μL</u>
Total volume	10 μ L

-Prepare a negative control reaction lacking *Bsa*I or diluted annealed oligo.

-Mix the reactions by pipetting. Incubate the reactions for 5-30 minutes at 37°C.

Note: Cleaning up the reactions with Zymo columns usually increases the number of colonies.

3. *E. COLI* TRANSFORMATION AND ANALYSIS OF TRANSFORMANTS

-Transform 1-4 μ l of the digestion-ligation reaction into an *E. coli* strain that doesn't have *ccd*B resistance (e.g. DH10B, TOP10, ...) to do counter-selection.

-Pick two colonies/construct, grow LB-Kan (100 mg/ml) cultures and purify plasmids.

-Sequence with appropriate primers: M13-F

(CCCAGTCACGACGTTGTAAAACGACGG) and M13-R

(CAGAGCTGCCAGGAAACAGCTATGACC) for *pENTR*-based vectors; attB1

(ACAAGTTTGTACAAAAAAGCAGGCT) and attB2

(ACCACTTTGTACAAGAAAGCTGGGT) primers for *pMDC32B*-, *pMDC123SB*- or *pFK210B*-based vectors).

Table I: *BsaI/ccdB*-based ('B/c') vectors for direct cloning of amiRNAs and syn-tasiRNAs.

Vector	Small RNA class	Bacterial antibiotic resistance	Plant antibiotic resistance	GATEWAY use	Backbone	Promoter	Terminator	Plant species tested
<i>pENTR-AtMIR390a-B/c</i>	amiRNA	Kanamycin	-	Donor	<i>pENTR</i>	-	-	-
<i>pFK210B-AtMIR390a-B/c</i>	amiRNA	Spectomycin	BASTA	-	<i>pGreen III</i>	<i>CaMV 35S</i>	<i>rbcS</i>	<i>A. thaliana</i>
<i>pMDC123SB-AtMIR390a-B/c</i>	amiRNA	Kanamycin	BASTA	-	<i>pMDC123</i>	<i>CaMV 2x35S</i>	-	<i>A. thaliana</i> <i>N. benthamiana</i>
<i>pMDC32B-AtMIR390a-B/c</i>	amiRNA	Kanamycin Hygromycin	Hygromycin	-	<i>pMDC32</i>	<i>CaMV 2x35S</i>	<i>nos</i>	<i>A. thaliana</i> <i>N. benthamiana</i>
<i>pENTR-AtTAS1c-B/c</i>	syn-tasiRNA	Kanamycin	-	Donor	<i>pENTR</i>	-	-	-
<i>pMDC123SB-AtTAS1c-B/c</i>	syn-tasiRNA	Kanamycin Hygromycin	BASTA	-	<i>pMDC123</i>	<i>CaMV 2x35S</i>	<i>nos</i>	<i>N. benthamiana</i>
<i>pMDC32B-AtTAS1c-B/c</i>	syn-tasiRNA	Kanamycin Hygromycin	Hygromycin	-	<i>pMDC32</i>	<i>CaMV 2x35S</i>	<i>nos</i>	<i>A. thaliana</i> <i>N. benthamiana</i>
<i>pENTR-AtTAS1c-D2-B/c</i>	syn-tasiRNA	Kanamycin	-	Donor	<i>pENTR</i>	-	-	-
<i>pMDC32B-AtTAS1c-D2-B/c</i>	syn-tasiRNA	Kanamycin Hygromycin	Hygromycin	-	<i>pMDC32</i>	<i>CaMV 2x35S</i>	<i>nos</i>	<i>N. benthamiana</i>